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

## COMFORT Epicondylitis-Spange

### COMFORT Epicondylitis Brace



**REF** 141 500 NET

ILLUSTRATIONEN  
FIGURES

1

2

SPRACHEN  
LANGUAGES

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	03
<b>EN</b>	English	Instructions for use .....	07

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Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### INTENDED PURPOSE

This medical device is a brace for relief of the muscle origins on the ulnar and/or radial epicondyle.

### INDICATIONS

Epicondylitis, chronic or recurring epicondylopathy, tendopathies, excessive strain, so-called tennis/golf arm.

### CONTRAINDICATIONS

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

### APPLICATION RISKS / IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

## **FITTING INSTRUCTIONS**

Put the brace on approx. 5 cm in front of the elbow joint. This position is normally equivalent to the highest circumference of the lower arm.

### **1 Epicondylitis radialis (tennis arm)**

In case of an epicondylitis radialis, the wide side of the brace should be placed on the upper muscle belly of the finger-stretching musculature with a pad. The studded side of the pad must be aligned in the direction of the wrist here.

### **2 Epicondylitis ulnaris (golfer arm)**

In case of an epicondylitis ulnaris, the wide side of the brace should be placed on the upper muscle belly of the finger flexor musculature with a pad. The studded side of the pad must be aligned in the direction of the wrist here.

Pull the restraints so tight that when strained – for example when closing the fist – pressure from the pad can be felt on the respective muscle belly. Ensure that no congestion occurs in the area which is located further away from the aid. If necessary, loosen the Velcro restraints.

## **REMOVAL**

To remove the epicondylitis brace, loosen the Velcro tape, increase the circumference and remove the brace in the direction of the hand.

## **MATERIAL COMPOSITION**

Polyamide (PA), polyurethane (PU), cotton (CO), polyester (PES), elasthane (EL)

The sewn in textile label provides the precise material composition.

## CLEANING INFORMATION

-  Hand-washing 30 °C
-  Do not bleach
-  Do not dry in a tumble dryer
-  Do not iron
-  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

## GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## USEFUL LIFE / LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link:  
[www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## DISPOSAL

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## DECLARATION OF CONFORMITY

The manufacturer confirms that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Status: 01.2025

 Medical device |  Single patient – multiple use

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PDF: [www.sanida.info](http://www.sanida.info)

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**NetworCare GmbH**

Goethestraße 5-7

09119 Chemnitz

Germany

Tel. 09421 8007 480

[www.sanida.info](http://www.sanida.info)

**BORT GmbH**  
Am Schweizerbach 1  
D-71384 Weinstadt  
[www.bort.com](http://www.bort.com)